INSTRUCTION FOR USE

A rapid test for detection of coronavirus SARS-CoV-2 antigen FOR PROFESSIONAL USE ONLY - NOT FOR HOME TESTING

REF: COVID19Speed-Antigen Test BSD_0503-25

INTENDED USE

COVID19Speed-Antigen Test is an immunochromatographic assay designed for the qualitative detection of specific SARS-CoV-2 antigen in human nasopharyngeal swab specimen. It is intended to be used for the detection of SARS-CoV-2 virus infection. The results obtained should not be the sole determinant for clinical decision.

CLINICAL INFORMATION

Coronaviruses are enveloped RNA viruses that are distributed broadly among humans, other mammals, and birds that cause respiratory, enteric, hepatic, and neurologic diseases (1, 2). Six coronavirus species are known to cause human disease (1, 2, 3). Four viruses — 229E, OC43, NL63, and HKU1 — are prevalent and typically cause common cold symptoms in immunocompetent individuals (1, 3). The two other strains — severe acute respiratory syndrome coronavirus (SARS-COV) and Middle East respiratory syndrome coronavirus (MERS-COV) — are zoonotic in origin and have been linked to sometimes fatal illness (4, 5). Coronaviruses are zoonotic, meaning they are transmitted between animals and people. It is spread by human-to-human transmission via droplets or direct contact, and infection has been estimated to have a mean incubation period of 6 days and a basic reproduction number of 2.2-3.6 (4, 6, 7). Common signs of infection include respiratory symptoms, fever, cough, shortness of breath and breathing difficulties. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure and even death (2, 3). A novel coronavirus (2019-nCoV, renamed as SARS-CoV-2 by WHO) was identified in early January 2020 as the infectious agent causing an outbreak of viral pneumonia in Wuhan, China, where the first cases had their symptom onset in December 2019 (8). It has since become an international health concern and declared a pandemic by WHO (9). The health emergency was declared to implement the resources necessary to control the infection (barrier measures, screening, treatment). In the absence of vaccine and specific treatment, the containment of the epidemics relies mainly on rapid identification and isolation of COVID-19 patients to prevent further spread of the virus (10, 11). Early diagnosis is important in this context, not only for the diagnosis and possible virological monitoring of hospitalized patients but also to ensure the health security of caregivers, first responders, and the general population (12). Therefore, development of point of care (POC) was recommended by WHO (8, 9, 13). Serological antibody-detection assays do not fulfil the requirement of the early detection (after infection) as the average incubation of 3-5 days is too short for development of an immune response. Real-time reverse transcriptase-polymerase chain reaction (RTPCR) is the main IVD assay used for COVID-19 and takes a few hours. Patients with Ct value >33 are no longer contagious (10, 11). The availability of a cost-effective, rapid POC antigen tests will play a critical role in the fight against COVID-19.

PRINCIPLE

COVID-19SPEEDAntigen Test is a colloidal gold enhanced double antibody sandwich immunoassay for qualitative determination of SARS-CoV-2 antigen in human nasal swabs. Antibodies to SARS-CoV-2 are immobilized in the test region (T) on nitrocellulose membrane. The target antigen in the specimen, if present, will react with gold conjugated anti-SARS-CoV-2 antibody and form an antibody-antigen complex. As this complex migrates along the length of the nitrocellulose membrane, it is captured by pre-coated anti-SARS-CoV-2 antibody located in the test (T) region on device test area (gold-SARS-CoV-2 Antibodies)-(SARS-CoV-2 Antigen)-(SARS-CoV-2 Antibodies), causing a pale to dark pink-purple band. If the specimen does not contain SARS-CoV-2 antigen, no colored line will appear in the test line region, indicating a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

MATERIAL PROVIDED

- 25 swabs
- 25 pcs of COVID-19 Rapid Antigen Test cassette packed in individually aluminum pouch with desiccant
- One bottle of sample Extraction Buffer
- 25 sets of Extraction tube and dripper
- Instruction for use (product insert)

MATERIALS REQUIRED BUT NOT PROVIDED

- Personal Protective Equipment, Timer, Biohazard container

STORAGE AND STABILITY

Store at 4-30°C, do not freeze. Keep the test device sealed until used. Keep away from direct sunlight, moisture and heat. When stored in a refrigerator, all kit components must be brought to room temperature (15-30 °C) for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.

LIMITATION OF THE TEST


- This product is designed for use with nasopharyngeal swab specimen only.
- This test detects the presence of specific SARS-CoV-2 antigen in the specimen and should not be used as the sole criterion for the diagnosis of a SARS-CoV-2 viral infection. The results obtained should be interpreted in conjunction with other diagnostic results and clinical information available to the physician.
- The test is a qualitative assay and it is not for quantitative determination of antigen concentration level. The intensity of the band does not have linear correlation with the antigen titer of the specimen.
- The performance of the test depends on antigen load in the specimen and may not correlate with PCR performed on the same specimen.
- Inadequate or inappropriate specimen collection, storage, and transportation may cause false negative results.
- If the test result is negative but symptoms persist, and a SARS-CoV-2 infection suspicion still exists, it is recommended to retest later or proceed follow-up testing using other clinical methods, for example RT-PCR method.
- A negative result at any time does not exclude the possibility of an early infection of SARS-CoV-2 virus.
- A negative result can occur if the quantity of SARS-CoV-2 virus antigen present in the specimen is below the detection limits of the assay.

WARNINGS AND PRECAUTIONS

- For professional in vitro diagnostics use only. Not for Home Testing.
- This product insert must be strictly followed in order to produce accurate test results.
- Keep the test device sealed until use. Perform the test within the 15 minutes after removing the test device from the foil pouch.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use the kit contents beyond the expiration date. The shelf life of the kit is as indicated on the outer package.
- The Buffer bottle may be opened and resealed for each assay. The Buffer cap should be firmly sealed between each use. The Buffer is stable until expiration date if kept at 2-30 °C.
- Avoid cross-contamination by using a new specimen extraction tube and dripper for each specimen preparation.
- Handle all specimens as being potentially infectious. Dispose all materials that come in contact with the specimen as infectious waste following the standard procedure.
- Wear protective clothing such as lab coat, disposable gloves and eye protection gear when specimens are being tested.
- Wipe any spills of sample promptly with 1% sodium hypochlorite solution or other disinfectant.
- Do not reuse test device. The provided sterilized swabs in the package should be used only for nasopharyngeal specimen collection.
- Direct swab specimens should be tested immediately after collection. If immediate testing is not possible, the swab specimen can be kept in an extraction tube filled with extraction buffer at room temperature (15-30 °C) for up to 30 minutes prior to testing.

TEST PROCEDURE (Refer to Figure)

SAMPLE COLLECTION


- Handle all specimens as being potentially infectious. Dispose of the used device according to your local regulations and biohazard waste disposal protocol. Dispose all materials that come in contact with the specimen as infectious waste.
- To collect nasopharyngeal swab specimen, it is important to obtain as much secretion as possible. Gently insert a sterile nasopharyngeal swab into the nasal cavity where the secretions are the most. Rotate the swab softly close to the inner wall of the nasal cavity 3 times and withdraw the swab. Make sure the tip of the swab is wet.


SAMPLE PREPARATION

- Add 200 µl (8 drops) of sample extraction buffer provided with the test kit to the extraction tube and then immerse the swab head in the tube completely. Mix the solution by rotating the swab vigorously and forcefully against the side of the tube at least 10 times (while submerged) and squeeze the tube 5 times by hand to ensure the sample on the swab is fully elutes into the sample extraction buffer.
- Squeeze the swab head along the inner wall of the extraction tube to keep the solution in the tube as much as possible. Discard the swab and cover the dripper head to mix the sample solution thoroughly.
- Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged periods. Samples are stable up to 30 mins when kept in the extraction buffer provided with the test kit.
- Samples should not be inactivated.
- If the sample has to be stored and transported, the swab should be stored in a clean, dry and closed container without media. The swab sample can be stored refrigerated at 2-8°C, or at room temperature (15-30°C) for up to 8 hours prior to testing. Dilution of the sample is not recommended, as dilution may cause decrease in test sensitivity.

SAMPLE PREPARATION


ASSAY PROCEDURE

INTERPRETATION OF RESULTS

Read the results within 5 - 15 minutes. Do not read the results after 20 minutes. Not following this procedure can lead to inaccurate results.

- **Negative:** Only control line (C) is visible. No SARS-CoV-2 virus antigen is detected. The result does not exclude COVID-19 infection (see the limitations section).
- **Positive:** Colored bands appear at the control line (C) and test line (T). The test is positive for SARS-CoV-2 virus antigen.
- **Invalid:** Control line (C) is absent. If this occurs, the assay should be repeated using a new test cassette.

PERFORMANCE CHARACTERISTICS

Preliminary Study - Performance data was calculated from a study of individuals suspected of exposure to COVID-19 with symptoms in the last 7 days, using a total of 39 known positive and 101 known negative samples; statistical analysis reported the comparison results are summarized in table below:

		RT-PCR		
		POS	NEG	TOTAL
BioSpeedia COVID-19Speed-Ag Test	POS	37	1	38
	NEG	2	100	102
	TOTAL	39	101	140
		Sensitivity: 0,95 CI (95%): 0,83-0,99	Specificity: 0,99 CI (95%): 0,95-0,99	Overall Agreement: 0.99

- Stratification of the positive specimens' post exposure between 0-3 days has a sensitivity of 100.0% (n=20) and 4-7 days has a sensitivity of 89.5% (n=19).
- Positive agreement of the BioSpeedia COVID-19Speed-Ag Test is higher with samples of CT values ≤32 with a sensitivity of 97.4%.

Analytical Specificity and Cross-reactivity - The BioSpeedia COVID-19Speed-Ag Test was evaluated with a total of 44 bacterial and viral isolates. Bacterial isolates were evaluated at a concentration between 10^7 and 10^9 org/mL. Viral isolates were evaluated at a concentration of at least 10^4 – 10^8 TCID₅₀/mL. Adenovirus 18 and Parainfluenza virus 3 were tested at 10^2 TCID₅₀/mL and 14 influenza viruses. There was no cross-reaction with the antibody/ diseases stated specimen tested.

Interfering Substances - Whole blood, and several over-the-counter (OTC) products and common chemicals were evaluated. No significant interference observed for the substances tested at the concentrations used [Whole blood (2%) - Three types of OTC mouthwashes (25%) - Three types of OTC throat drops (25%) - Three types of OTC nasal sprays (10%) - 4-Acetamidophenol (10 mg/mL) - Acetylsalicylic Acid (20 mg/mL) - Chlorpheniramine (5 mg/mL) - Dextromethorphan (10 mg/mL) - Diphenhydramine (5 mg/mL) - Ephedrine (20 mg/mL) - Guaiaicol glyceryl ether (20 mg/mL) - Oxymetazoline (10 mg/mL) - Phenylephrine (100 mg/mL) - Phenylpropanolamine (20 mg/mL)].

Hook effect - There is no hook effect at $1.0 \times 10^{5.9}$ TCID₅₀/ml of SARS- CoV-2.

Temperature Stability Study - Stable at 45°C for 56 days, shelf life at least 18 months from manufacturing date under 4-30°C is justified.

Shipping Stability Study - Stable when shipping in its original packing condition, at ambient temperature condition.

Repeatability & Reproducibility - Repeatability & Reproducibility was established using in-house reference panels containing negative specimens and a range of positive specimens. There were no differences observed within-run, between-run, between-lots, between-sites, and between-days.

WARRANTY AND LIMITED LIABILITY

Failure to follow the assay procedure may derive inaccurate results. In such event, the manufacturer disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and the fitness for use. The manufacturer will not be liable for any damage caused by misuse, improper handling and storage, non-compliance with warnings and procedures, damage caused by events occurring after the product is released, failure to ensure the product is in proper condition before use, or any warranty given by independent distributor.

REFERENCES

1. Weiss SR, Leibowitz JL. Coronavirus pathogenesis. *Advances in Virus Research* (2011), 81:85-164.
2. Rothan HA, Byrareddy SN. The epidemiology and pathogenesis of coronavirus disease (COVID-19) outbreak. *J Autoimmun.* 2020; Feb 26:102433. doi: 10.1016/j.jaut.2020.102433.
3. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. *Trends in Microbiology* (2016), 24:490-502.
4. Guo YR, Cao QD, Hong ZS, et al. The origin, transmission and clinical therapies on coronavirus disease 2019 (COVID-19) outbreak-an update on the status. *Mil Med Res.* 2020; Mar 13; 7(1):11.doi:10.1186/s40779-020-00240-0.
5. Lai CC, Shih TP, Ko WC, et al. Severe acute respiratory syndrome coronavirus 2 (SARS-Cov-2) and coronavirus disease-2019 (COVID-19): The epidemic and the challenges. *Int J Antimicrob Agents.* 2020; Mar 55(3): 105924.doi: 10.1016/j.ijantimicag.2020.105924.
6. CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
7. CDC. Duration of Isolation and Precautions for Adults with COVID-19. (2020).
8. World Health Organization (WHO). WHO Statement Regarding Cluster of Pneumonia Cases in Wuhan, China. Beijing: WHO; 9 Jan 2020. [Accessed 26 Jan 2020]. <https://www.who.int/china/news/detail/09-01-2020-who-statement-regarding-cluster-of-pneumonia-cases-in-wuhan-china>.
9. World Health Organization (WHO). Coronavirus. <https://www.who.int/health-topics/coronavirus>.
10. CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
11. CDC. Duration of Isolation and Precautions for Adults with COVID-19. (2020).

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RÉFÉRENCES

- Weiss SR, Leibowitz JL. Coronavirus pathogenesis. *Advances in Virus Research* (2011), 81:85-164.
- Rothan HA, Byrareddy SN. The epidemiology and pathogenesis of coronavirus disease (COVID-19) outbreak. *J Autoimmun.* 2020; Feb 26:102433. doi: 10.1016/j.jaut.2020.102433.
- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. *Trends in Microbiology* (2016), 24:490-502.
- Guo YR, Cao QD, Hong ZS, et al. The origin, transmission and clinical therapies on coronavirus disease 2019 (COVID-19) outbreak-an update on the status. *Mil Med Res.* 2020; Mar 13; 7(1):11.doi:10.1186/s40779-020-00240-0.
- Lai CC, Shih TP, Ko WC, et al. Severe acute respiratory syndrome coronavirus 2 (SARS-Cov-2) and coronavirus disease-2019 (COVID-19): The epidemic and the challenges. *Int J Antimicrob Agents.* 2020; Mar 55(3): 105924.doi: 10.1016/j.ijantimicag.2020.105924.
- CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
- CDC. Duration of Isolation and Precautions for Adults with COVID-19. (2020).
- World Health Organization (WHO). WHO Statement Regarding Cluster of Pneumonia Cases in Wuhan, China. Beijing: WHO; 9 Jan 2020. [Accessed 26 Jan 2020]. <https://www.who.int/china/news/detail/09-01-2020-who-statement-regarding-cluster-of-pneumonia-cases-in-wuhan-china>.
- World Health Organization (WHO). Coronavirus. <https://www.who.int/health-topics/coronavirus>.
- CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
- CDC. Duration of Isolation and Precautions for Adults with COVID-19. (2020).
- In Vitro Diagnostic Assays for COVID-19: Recent Advances and Emerging Trends (Sandeep Kumar Vashist, 2020 April 05: diagnostics).

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
REFERENCIAS

25. Weiss SR, Leibowitz JL. Coronavirus pathogenesis. *Advances in Virus Research* (2011), 81:85-164.
26. Rothan HA, Byrareddy SN. The epidemiology and pathogenesis of coronavirus disease (COVID-19) outbreak. *J Autoimmun.* 2020; Feb 26:102433. doi: 10.1016/j.jaut.2020.102433.
27. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. *Trends in Microbiology* (2016), 24:490-502.
28. Guo YR, Cao QD, Hong ZS, et al. The origin, transmission and clinical therapies on coronavirus disease 2019 (COVID-19) outbreak-an update on the status. *Mil Med Res.* 2020; Mar 13; 7(1):11.doi:10.1186/s40779-020-00240-0.
29. Lai CC, Shih TP, Ko WC, et al. Severe acute respiratory syndrome coronavirus 2 (SARS-Cov-2) and coronavirus disease-2019 (COVID-19): The epidemic and the challenges. *Int J Antimicrob Agents.* 2020; Mar 55(3): 105924.doi: 10.1016/j.ijantimicag.2020.105924.
30. CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
31. CDC. Duration of Isolation and Precautions for Adults with COVID-19. (2020).
32. World Health Organization (WHO). WHO Statement Regarding Cluster of Pneumonia Cases in Wuhan, China. Beijing: WHO; 9 Jan 2020. [Accessed 26 Jan 2020]. <https://www.who.int/china/news/detail/09-01-2020-who-statement-regarding-cluster-of-pneumonia-cases-in-wuhan-china>.
33. World Health Organization (WHO). Coronavirus. <https://www.who.int/health-topics/coronavirus>.
34. CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
35. CDC. Duration of Isolation and Precautions for Adults with COVID-19. (2020).
36. In Vitro Diagnostic Assays for COVID-19: Recent Advances and Emerging Trends (Sandeep Kumar Vashist, 2020 April 05: diagnostics).

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BioSpeedia International S.A.R.L