



A spin-off from Institut Pasteur

## END-TO-END SOLUTIONS for the screening of SARS-CoV-2 (COVID-19)

BioSpeedia is a spin-off from Institut Pasteur, the premier non-profit foundation that was founded in 1887 and has since originated 10 Nobel Prizes.

Launched in 2011, the company specializes in rapid tests across various infectious diseases.


### BioSpeedia Overview



BioSpeedia is a French company headquartered at the innovation campus at St Etienne and has privileged access to **CHU of St Etienne, GIMAP and Ecole des Mines** and has a license from the Institut Pasteur granting access to major reagents.



BioSpeedia provides **both Antigen and Antibody** rapid testing kits that are **fast, accurate, and totally secured**.



BioSpeedia provides **an extraction kit as well as a Multiplex, RT –qPCR Kit**.



BioSpeedia can provide up to **5 million test kits per month for each test** through our existing manufacturing facilities.





**end-to-end solutions**  
for COVID-19 screening

**Rapid tests for  
COVID-19**



**Antibody tests**

**Why are antibody  
tests important**



**How to use the tests**


## BioSpeedia Product Covid-19 Portfolio



**Rapid antigen test** – 5 to 15 minutes – for the qualitative detection of specific SARS-CoV-2 antigen in human nasopharyngeal swab specimen.  
**Saliva and Sputum antigen test under development.**



**Rapid antibody test** – 5 to 10 minutes – for the qualitative detection of specific IgM and IgG antibodies against SARS-CoV-2, in human whole blood, serum or plasma.



**Real-time PCR testing kit** for the detection of 2019-nCoV genome; in this one-step real-time RT-PCR, reverse transcription of viral RNA is combined with the qPCR step in a single tube reaction.

## Difference between an antibody test and an antigen test



### Antigen test

Antigen tests use nasal or throat swabs to detect **active COVID-19 infection**.

Antigen tests detect specific proteins on the surface of the virus. Antigen tests are similar to PCR tests but are **quicker and less expensive** than a PCR test with also a **high sensitivity/specificity**.



### Antibody test

Serology (blood-based) antibody tests detect presence of antibodies for **an active or past COVID-19 infection**.

Antibody tests show that you **might still be infected (IgM) or were previously infected (IgG) by the virus**, even if you never showed any signs or symptoms (asymptomatic).



## BioSpeedia Rapid Antigen Test

IVD

**Description:** This is an immunochromatographic assay designed for the qualitative detection of specific SARS-CoV-2 antigen in human nasopharyngeal swab specimen. It is intended to be used by the professionals as a screening test and as an aid in the early diagnosis of SARS-CoV-2 infection.

### Features

Platform	Immunochromatographic Assay
Format	Cassette
Detection	Specific SAR-CoV-2 antigen
Specimen	Nasopharyngeal swabs
Assay time	5-15 minutes
Shelf life	18 months

**Clinical performance:** A study of 324 patients shows a sensitivity of 97.4%\* and a specificity of 99%. The study was done in 2 French hospitals and also shows higher cut-off versus similar tests.

Results	PCR Positive	Test Negative	Total	Sensitivity (PPA, Positive Percent Agreement)	Specificity (NPA, Negative Percent Agreement)
Positive	107	3	110	94.7% (107/113)* * Sensitivity = 97.4% for CT < 32	99% (208/211)
Negative	6	208	214		
Total	113	211	324		

### Procedures (Simplified Version)

1



Add 200 µl (~ 8 drops) extraction buffer to the tube

2



Immerse the swab head and mix vigorously at least 10 times

3




Squeeze liquid from swab and discard the swab head

4



Cover the dropper head and mix thoroughly


5



Add 50 µL (2 drops) to the sample well. Read result in 10-15 minutes

### Result interpretation

C : Control Line  
T : Test Line





**Definition:** This is an immunochromatographic assay designed for the qualitative detection of specific IgM and IgG antibodies against SARS-CoV-2, in human whole blood, serum or plasma. It is intended to be used by the professionals as a screening test and as an aid in the early diagnosis of SARS-CoV-2 infection.

## Features

Platform	Immunochromatographic Assay
Format	Cassette
Detection	Specific IgM and IgG antibodies against SARS-CoV-2
Specimen	Whole blood, serum or plasma
Specimen volume	20 µl
Assay time	5-10 minutes
Shelf life	18 months

**Clinical performance:** COVID-19 Rapid IgG/IgM Test has been validated by 5 hospitals located in 3 European countries (Italy, France, Germany). A total of 564 patients, all confirmed positive for SARS-CoV-2 RNA by RT-PCR, including 549 symptomatic and 15 asymptomatic patients, were included for the sensitivity analysis. Results below:

Results (IgM+IgG)	PCR Positive	Test Negative	Total	Sensitivity (PPA, Positive Percent Agreement)	Specificity (NPA, Negative Percent Agreement)
Positive	392	4	396	92.2% (392/425)	98.1% (211/215)
Negative	33	211	244		
Total	425	215	640		


## Procedures (Simplified Version)

1




Add 20 µL of whole blood, serum or plasma to the sample well. Make sure there is no air bubble.

2



Add 1 drop of buffer to the same sample well. Read result within 15 minutes.

## Result Interpretation





## Why use antibody testing?



Conduct epidemiological studies; plan for social distancing



Find good candidates for convalescent plasma therapy



Identify children with MIS



Reassess medication for immune compromised patients



Identify *first-in-line* candidates for vaccination



Identify high-risk patients among asymptomatic populations

## Why antibody testing is important



In communities with lower prevalence (5 – 25%), infection risk is very high. **Antibody tests** allow governments and organizations to remain alert on that potential risk and adjust measures and resources to **prevent** significant **loss of lives and businesses**.



Around 25% (approx. 1.8 billion people) of the world's population has **underlying medical conditions** including old-age, obesity, severe cardiac, liver, lung, asthma, diabetes, kidney. Antibody tests allow **high-risk asymptomatic** patients to know if they've been previously infected by COVID-19 and whether their underlying conditions have been adversely impacted and/or if any change is required in their prescription medications.



Periodic anti-body testing allows organizations to address their exposure to COVID-19 and assess whether additional distancing measures should be taken.



# Combo 2019-nCoV Viral Extraction and Real-time RT-qPCR Kit

## A combination of:

- GF-1 Viral Nucleic Acid Extraction Kit and 2019-nCoV Multiplex RT-qPCR Kit
- A-Z COVID-19 testing solution from COVID-19 RNA extraction till RT-qPCR

**Easy and convenient:** Able to complete the entire process from nucleic acid purification to RT-qPCR assay within 3.5 hours

## Nucleic Acid Extraction

### GF-1 Viral Nucleic Acid Extraction Kit (Mini-Prep Kit)

- Rapid and efficient purification of viral DNA/RNA from samples such as serum, plasma, body fluid, supernatant of virus-infected cell culture or supernatant of virus-infected swabs and/or nasopharyngeal swab
- Efficient recovery of highly pure DNA or RNA
- Purified DNA/RNA is suitable for downstream reactions, such as RT-qPCR and library prep for next generation sequencing
- Eliminates the need of toxic or organic-based extraction

## One Step RT-qPCR Master Mix

### One Step AtTaq RT-qPCR Master Mix

- Formulation optimized with Hot Start Taq, M-MULV enzyme and enhancers
- ROX dye is supplied separately
- Compatible with most of the real-time PCR platform
- Multiplex up to 3 targets

Description	Pack Size
Combo 2019-nCoV Viral Extraction and Real-time RT-qPCR Kit	100 preps
GF-1 Viral Nucleic Acid Extraction Kit (Proteinase K & Carrier RNA included)	100 preps
2019-nCoV Multiplex RT-qPCR Kit	100 preps
One Step AtTaq RT-qPCR Master Mix	100 preps

COVID 19






# Combo 2019-nCoV Viral Extraction and Real-time RT-qPCR Kit



## Multiplex RT-qPCR Assay

### 2019-nCoV Multiplex RT-qPCR Kit

- Multiplex detection in a single tube: N gene, RdRp gene and E gene
- Ready-to-use solution with all necessary components (primer, probe and master mix) in a single tube
- Internal control (IC): human gene is included
- Minimize technical error
- Enable immediate reaction set up with short processing time



Well/Sample	Ct Value			
	FAM	CY5	Texas Red	Hex
Pos Control 1	14.36	12.43	12.87	15.01
Pos Control 2	13.94	12.39	12.45	14.07
Pos Control 3	14.43	12.68	13.42	15.11
Pos Control 4	14.24	12.5	13.36	14.56
Pos Control 5	14.41	12.1	15.01	14.11
Pos Control 6	13.62	12.01	14.06	14.18
NTC	No Ct	No Ct	No Ct	No Ct

The graphs show the consistency of the kit performance for different genes in 2019-nCoV Multiplex RT-qPCR Kit. The kit was tested for stability from week one until week six for ready-to-use One Step RT-qPCR pre-mixed with four different genes of primers and probes in the real-time master mixes. The graphs and Ct values show that the kit has stable and consistent performance.



## A Biospeedia SICPA® Solution



Biospeedia is partnering with SICPA® in developing a **secured and certified health pass** featuring COVID-tests results to allow secured traveling.

**SICPA®** is a trusted global provider of security inks as well as a leader in secured identification, traceability and authentication solutions.




**CERTUS™ myHealth Pass** provides guaranteed and secured proof of test status in both digital and material form. It protects against fraud thanks to a cutting edge blockchain-secured digital seal, based on the Guardtime KSI® Blockchain.



**Biospeedia** Test status data is revocable and thus always up-to-date. It respects personal privacy concerns while providing incontrovertible assurances of health status to the aviation and border authorities.

The credential is authenticated by scanning a tamper-proof QR-code with a smartphone or a computer. The verification is independent of the document issuer. No personal or sensitive information is stored, either directly or encrypted, in any data base or on the blockchain and no extra-infrastructure is required.

A Biospeedia SICPA® Solution



Passengers can take a test prior to the flight (**72-24 hours**) at a clinic or hospital by healthcare professionals (doctors, pharmacists, certified nurses)



Results in 15 minutes are saved in "myHealth Pass"



Passengers can take a test prior to the flight (**<4 hours**) at the airport done by onsite certified professionals



Results in 15 minutes are saved in "myHealth Pass"



Positive Result: quarantine and re-test. Can't board flight



Negative Result: authorities can access "myHealth Pass" and allow passengers to board the flight



Positive Result: quarantine and re-test. Can't board flight



## ORDERING INFORMATION



### BioSpeedia Rapid Antigen Test



Cat. No	Product Description	Test Specimen	Product Format	Intended Use	Storage Temp	Packing Size
BSD_0503-25	COVID19Speed-Antigen (with Nasopharyngeal Swab)	Nasopharyngeal swab sample	Cassette	Coronavirus disease	4-30°C	25 T/ Kit
BSD_0503-10	COVID19Speed-Antigen (with Nasopharyngeal Swab)	Nasopharyngeal swab sample	Cassette	Coronavirus disease	4-30°C	10 T/ Kit



### BioSpeedia Rapid Antibody Test



Cat. No	Product Description	Test Specimen	Product Format	Intended Use	Storage Temp	Packing Size
BSD_0501	COVID19SEROspeed-IgM-IgG (with disposable dropper, alcohol swab, lancet)	Whole blood, serum or plasma	Cassette	Coronavirus disease	4-30°C	25 T/ Kit



### BioSpeedia Viral Nucleic Acid Extraction Kit



Cat. No	Product Description	Test Specimen	Intended Use	Storage Temp	Packing Size
BSD_0510	Viral Nucleic Acid Extraction Kit	Serum, plasma, body fluid, virus-infected cell culture supernatant, VTM buffers consists of virus infected swabs	Isolation and purification of Viral nucleic Acid (DNA/RNA)	4-30°C	100 Preps



### BioSpeedia 2019-nCoV Multiplex RT-qPCR Kit



Cat. No	Product Description	Test Specimen	Intended Use	Storage Temp	Packing Size
BSD_0511	2019-nCoV Multiplex RT-qPCR Kit	Extracted RNA from nasopharyngeal & oropharyngeal swab	Detection of SARS-CoV-2's RNA from nasopharyngeal and oropharyngeal swabs	-20°C	100 Preps

